

Product datasheet for **RC203892A1V**

Human **GNGT2 (NM_031498) AAV Particle**

Product data:

Product Type: AAV Particles
Product Name: Human GNGT2 (NM_031498) AAV Particle
Tag: Myc-DDK
Symbol: GNGT2
Synonyms: G-GAMMA-8; G-GAMMA-C; GNG9; GNGT8
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203892 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCC**CCAGGATCTCAGCGAGAAGGACCTGTTGAAGATGGAGGTGGAGCAGCTGAAGAAAGAAGTGAAAA**
ACACAAGAATTCGATTTCAAAGCGGGAAAGGAAATCAAGGAGTACGTGGAGGCCCAAGCAGGAAACGA
TCCTTTTCTCAAAGGCATCCCTGAGGACAAGAATCCCTTCAAGGAGAAAGGTGGCTGTCTGATAAGC

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203892 protein sequence
 Red=Cloning site Green=Tags(s)

MAQDLSEKDLLKMEVEQLKKEVKNTRIPISKAGKEIKEYVEAQAGNDPFLKGIPEDKNPFKEKGGCLIS

TRTRPLE**QKLISEEDLAANDILDYKDDDDKV**

Species: Human
Serotype: AAV-2
ACCN: NM_031498
ORF Size: 207 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_031498.1](#), [NP_113686.1](#)

RefSeq Size: 1057 bp

RefSeq ORF: 210 bp

Locus ID: 2793

UniProt ID: [O14610](#)

Cytogenetics: 17q21.32

MW: 7.7 kDa