

Product datasheet for **RC203896L1V**

PCSK7 (NM_004716) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PCSK7 (NM_004716) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PCSK7 |
| Synonyms: | LPC; PC7; PC8; SPC7 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_004716 |
| ORF Size: | 2355 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC203896). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_004716.2 |
| RefSeq Size: | 3940 bp |
| RefSeq ORF: | 2358 bp |
| Locus ID: | 9159 |
| UniProt ID: | Q16549 |
| Cytogenetics: | 11q23.3 |
| Domains: | Peptidase_S8, P_protein |
| Protein Families: | Druggable Genome, Protease, Transmembrane |



[View online »](#)

MW: 86.2 kDa

Gene Summary: This gene encodes a member of the subtilisin-like proprotein convertase family, which includes proteases that process protein and peptide precursors trafficking through regulated or constitutive branches of the secretory pathway. It encodes a type 1 membrane bound protease that is expressed in many tissues, including neuroendocrine, liver, gut, and brain. The encoded protein undergoes an initial autocatalytic processing event in the ER and then sorts to the trans-Golgi network through endosomes where a second autocatalytic event takes place and the catalytic activity is acquired. This gene encodes one of the seven basic amino acid-specific members which cleave their substrates at single or paired basic residues. It can process proalbumin and is thought to be responsible for the activation of HIV envelope glycoproteins gp160 and gp140. This gene has been implicated in the transcriptional regulation of housekeeping genes and plays a role in the regulation of iron metabolism. A t(11;14)(q23;q32) chromosome translocation associated with B-cell lymphoma occurs between this gene and its inverted counterpart. [provided by RefSeq, Feb 2014]