

Product datasheet for **RC203901L4V**

CNRIP1 (NM_015463) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CNRIP1 (NM_015463) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CNRIP1
Synonyms:	C2orf32; CRIP-1; CRIP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015463
ORF Size:	492 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203901).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015463.1
RefSeq Size:	2037 bp
RefSeq ORF:	495 bp
Locus ID:	25927
UniProt ID:	Q96F85
Cytogenetics:	2p14
MW:	18.6 kDa


[View online »](#)

Gene Summary:

This gene encodes a protein that interacts with the C-terminal tail of cannabinoid receptor 1. Two transcript variants encoding different isoforms have been described for this gene.
[provided by RefSeq, Jul 2013]