

## Product datasheet for **RC203906L4V**

### **RBCK1 (NM\_031229) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RBCK1 (NM_031229) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBCK1
Synonyms:	C20orf18; HOIL-1; HOIL1; PBMEI; PGBM1; RBCK2; RNF54; UBCE7IP3; XAP3; XAP4; ZRANB4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_031229
ORF Size:	1500 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203906).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031229.2</a> , <a href="#">NP_112506.1</a>
RefSeq Size:	2778 bp
RefSeq ORF:	1533 bp
Locus ID:	10616
UniProt ID:	<a href="#">Q9BYM8</a>
Cytogenetics:	20p13
Domains:	RING, zf-RanBP
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 56.4 kDa

**Gene Summary:** The protein encoded by this gene is similar to mouse UIP28/UbcM4 interacting protein. Alternative splicing has been observed at this locus, resulting in distinct isoforms. [provided by RefSeq, Jul 2008]