

Product datasheet for **RC203912A1V**

Human TNFAIP8L1 (NM_152362) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TNFAIP8L1 (NM_152362) AAV Particle
Tag: Myc-DDK
Symbol: TNFAIP8L1
Synonyms: TIPE1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203912 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGACACCTTCAGCACCAAGAGCCTGGCTCTGCAGGCGCAGAAGAAGCTCCTGAGTAAGATGGCGTCCA
AGGCAGTGGTGGCCGTGCTGGTGGATGACACCAGCAGTGAGGTGCTGGATGAGCTGTACCGCGCCACCAG
GGAGTTACGCGCAGCCGCAAGGAGGCCGAGAAGATGCTCAAGAACCTGGTCAAGGTGGCCCTGAAGCTG
GACTGCTGCTGCGTGGGGACCAGCTGGGCGGTGAGGAGCTGGCGCTGCTGCGGCGCTTCCGCCACCGGG
CGCGCTGCCTGGCCATGACGGCCGTGAGTTCACCCAGGTGGACTTCACCTTCGACCGGCGCGTGTGGC
CGTCGGGCTGCTCGAGTGCCGCGACCTGCTGCACAGGCGTGGGTCCCACTGACCGCCAAGTCCAC
GGCCGCATCAACCACGTGTTGCGCCACCTAGCCGACTGCGACTTCCTGGCTGCGCTCTACGGCCCCGCGC
AGCCCTACCGCTCCCACTGCGCAGGATCTGCGAGGCGCTGGCCGGATGCTGGACGAGGGCAGCCTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203912 protein sequence
Red=Cloning site Green=Tags(s)

MDTFSTKSLALQAQKLLSKMASKAVVAVLVDDTSSEVLDEL YRATREFTRSRKEAQKMLKNLVKVALKL
GLLLRGDQLGGEELALLRRFRHRARCLAMTAVSFHQVDFTFDRRVLAVGLLECRDLLHQAVGPHLTAKSH
GRINHVFGHLADCDFLAALYGPAEYPYRSHLRRIEGLGRMLDEGSL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_152362
ORF Size:	558 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152362.1</u> , <u>NP_689575.1</u>
RefSeq Size:	3817 bp
RefSeq ORF:	561 bp
Locus ID:	126282
UniProt ID:	<u>Q8WVP5</u>
Cytogenetics:	19p13.3
MW:	20.9 kDa