

Product datasheet for **RC203912L3V**

TNFAIP8L1 (NM_152362) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TNFAIP8L1 (NM_152362) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TNFAIP8L1
Synonyms:	TIPE1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_152362
ORF Size:	558 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203912).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_152362.1 , NP_689575.1
RefSeq Size:	3817 bp
RefSeq ORF:	561 bp
Locus ID:	126282
UniProt ID:	Q8WVP5
Cytogenetics:	19p13.3
MW:	20.9 kDa
Gene Summary:	Acts as a negative regulator of mTOR activity.[UniProtKB/Swiss-Prot Function]


[View online »](#)