

## Product datasheet for **RC203923L2V**

### **RHBDF2 (NM\_024599) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RHBDF2 (NM_024599) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RHBDF2
Synonyms:	iRhom2; RHBDL5; RHBDL6; TEC; TOC; TOCG
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_024599
ORF Size:	2568 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203923).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024599.2</a> , <a href="#">NP_078875.2</a>
RefSeq Size:	3615 bp
RefSeq ORF:	2571 bp
Locus ID:	79651
UniProt ID:	<a href="#">Q6PJF5</a>
Cytogenetics:	17q25.1
Protein Families:	Protease, Transmembrane
MW:	96.7 kDa



[View online »](#)

**Gene Summary:**

Regulates ADAM17 protease, a sheddase of the epidermal growth factor (EGF) receptor ligands and TNF, thereby plays a role in sleep, cell survival, proliferation, migration and inflammation. Does not exhibit any protease activity on its own.[UniProtKB/Swiss-Prot Function]