

Human COUP TF1 (NR2F1) (NM_005654) AAV Particle

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	COUP TFI
Synonyms:	BBOAS; BBSOAS; COUP-TFI; COUPTFI; EAR-3; EAR3; ERBAL3; SVP44; TCFCOUP1; TFCOUP1
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203924 representing NM_005654 Red=Cloning site Blue=ORF Green=Tags(s)

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC203924 representing NM_005654 Red=Cloning site Green=Tags(s) MAMVVSSWRDPQDDVAGGNPGGPNPAAQAARGGGGAGEQQQQAGSGAPHTPQTGQPGAPATPGTAGDK GQGPPGSGQSQQHIECVVCGDKSSGKHYGQFTCEGCKSFFKRSVRRNLTYTCRANRNCPIDQHHRNQCQY CRLKKCLKVGMREAVQGRMPPTQPNPGQYALTNGDPLNGHCYLSGYISLLRAEPYPTSRYGSCMQP NNIMGIENICELAAARLLFAVEWARNIPFFDLQITDQVSLRLTWSELFVLNAAQCSMPLHVAPLLAAA GLHASPMSADRVVAFMDHIRIFQEVEKLKALHVDSAEYSCLKAIVLFTSDACGLSDAAHIESLQEKSSQC ALEEYVRSQYPNQPSRFGKLLRLPSLRTVSSSVIEQLFFVRLVGKTPJETLIRDMLLSGSSFNWPYMSI QCS TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005654
ORF Size:	1269 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005654.4</u>
RefSeq Size:	3210 bp
RefSeq ORF:	1272 bp
Locus ID:	7025
UniProt ID:	<u>P10589</u>
Cytogenetics:	5q15
MW:	46.6 kDa