

Product datasheet for **RC203924L4V**

COUP TF1 (NR2F1) (NM_005654) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	COUP TF1 (NR2F1) (NM_005654) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NR2F1
Synonyms:	BBOAS; BBSOAS; COUP-TFI; COUPTF1; EAR-3; EAR3; ERBAL3; SVP44; TCFCOUP1; TFCOUP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005654
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203924).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005654.4
RefSeq Size:	3210 bp
RefSeq ORF:	1272 bp
Locus ID:	7025
UniProt ID:	P10589
Cytogenetics:	5q15
Domains:	HOLI, zf-C4
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors


[View online »](#)

MW: 46.6 kDa

Gene Summary: The protein encoded by this gene is a nuclear hormone receptor and transcriptional regulator. The encoded protein acts as a homodimer and binds to 5'-AGGTCA-3' repeats. Defects in this gene are a cause of Bosch-Boonstra optic atrophy syndrome (BBOAS). [provided by RefSeq, Apr 2014]