

## Product datasheet for RC203925L4

### ZUFSP (ZUP1) (NM\_145062) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZUFSP (ZUP1) (NM_145062) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	ZUP1
Synonyms:	C6orf113; DUB; ZUFSP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203925).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_145062
ORF Size:	1734 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_145062.1</a>
<b>RefSeq Size:</b>	2226 bp
<b>RefSeq ORF:</b>	1737 bp
<b>Locus ID:</b>	221302
<b>UniProt ID:</b>	<a href="#">Q96AP4</a>
<b>Cytogenetics:</b>	6q22.1
<b>Domains:</b>	zf-C2H2
<b>MW:</b>	66 kDa
<b>Gene Summary:</b>	This gene encodes a protein containing zinc finger motifs and a cysteine peptidase domain. The encoded protein functions as a K63-specific de-ubiquitinating enzyme that specifically cleaves long K63-linked polyubiquitin chains in the middle of a chain (i.e. "endo cleavage") rather than by removing the terminal ubiquitin from a chain. This enzyme is thought to be involved in the regulation of DNA repair by cleaving K63-linked ubiquitin chains at repair foci. This protein is related to proteases for the ubiquitin-like modifiers Ufm1 (ubiquitin fold modifier 1) and Atg8/Gabarapl2, but does not have any activity on these modifiers. [provided by RefSeq, Mar 2018]

Product images:



Circular map for RC203925L4