

Product datasheet for **RC203942A1V**

Human NCE2 (UBE2F) (NM_080678) AAV Particle

Product data:

Product Type: AAV Particles
 Product Name: Human NCE2 (UBE2F) (NM_080678) AAV Particle
 Tag: Myc-DDK
 Symbol: NCE2
 Synonyms: NCE2
 Mammalian Cell Selection: None
 Vector: pAAV-AC-Myc-DDK (PS100089)
 ORF Nucleotide Sequence: >RC203942 representing NM_080678
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTAACGCTAGCAAGTAACTGAAGCGTGACGATGGTCTCAAAGGGTCCCGGACGGCAGCCACAGCGT
 CCGACTCGACTCGGAGGGTTTCTGTGAGAGACAAATTGCTTGTAAAGAGGTTGCAGAACTGAAGCTAA
 TTTACCTTGATACATGTAAAGTGCATTTTCCTGATCCAAACAAGCTTCATTGTTTTAGCTAACAGTAACC
 CCAGATGAGGGTTACTACCGGGTGGAAATTTTCAGTTTGAACTGAAGTCCCGATGCGTACAACATGG
 TGCTTCCCAAAGTGAAATGCCTGACCAAGATCTGGCACCCCAACATCACAGAGACAGGGGAAATATGCT
 GAGTTTATTGAGAGAACATTCAATTGATGGCACTGGCTGGGCTCCACAAGAACATTAAGGATGTCGTT
 TGGGGATTAACTCTTTGTTTACTGATCTTTGAATTTGATGATCCACTGAATATTGAAGCTGCAGAAC
 ATCATTGCGGGACAAGGAGGACTTCCGGAATAAAGTGGATGACTACATCAAACGTTATGCCAGA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203942 representing NM_080678
 Red=Cloning site Green=Tags(s)

MLTLASKLKRDDGLKGSRTAATASDSTRVSVRDKLLVKEVALEANLPCTCKVHFDPNKLHCFQLTVT
 PDEGYQGKQFETEVDPAYNMVPPKVKCLTKIWHPNITETGEICLSLLREHSIDGTGWAPRTLKDVV
 WGLNSLFTDLLNFDDPLNIEAAEHHLRDKEDFRNKVDDYIKRYAR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
 Serotype: AAV-2



[View online »](#)

ACCN:	NM_080678
ORF Size:	555 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_080678.1</u>
RefSeq Size:	1366 bp
RefSeq ORF:	558 bp
Locus ID:	140739
UniProt ID:	<u>Q969M7</u>
Cytogenetics:	2q37.3
MW:	20.9 kDa