

## Product datasheet for **RC203944A1V**

### Human MAX (NM\_145114) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human MAX (NM\_145114) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** MAX  
**Synonyms:** bHLHd4  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC203944 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGCGATAACGATGACATCGAGGTGGAGAGCGACGAAGAGCAACCGAGGTTTCAATCTGCGGCTGACA  
AACGGGCTCATCATAATGCACTGGAACGAAAACGTAGGGACCACATCAAAGACAGCTTTCACAGTTTGCG  
GGACTCAGTCCCATCACTCCAAGGAGAGAAGCTCTATTTCTCTTTTGGAAATTGTACTCCTGTCTT  
CATCGTCAAAGTTTGATGCAGAAATGCCACACCTTCATTTCAAGCTACCAAGTGCACAAGAAAAAAGAAT  
GCAAGATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC203944 protein sequence  
Red=Cloning site Green=Tags(s)

MSDNDDIEVESDEEQPRFQSAADKRAHHNALERKRRDHKDSFHSRLRDSVPSLQGEKLYFLFWKLCTPVL  
HRQSLMQKCHTFISSYQVHKKECKI

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_145114  
**ORF Size:** 288 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_145114.1</u>
<b>RefSeq Size:</b>	937 bp
<b>RefSeq ORF:</b>	291 bp
<b>Locus ID:</b>	4149
<b>UniProt ID:</b>	<u>P61244</u>
<b>Cytogenetics:</b>	14q23.3
<b>MW:</b>	11.5 kDa