

Product datasheet for **RC203963**

WDR61 (NM_025234) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR61 (NM_025234) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	WDR61
Synonyms:	REC14; SKI8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC203963 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCAACCAAGTACGGTATTCTCTTCAAACAAGAGCAAGCCCATGATGATGCCATTTGGTCAGTTGCTT
GGGGGACAAACAAGAAGGAAAACCTGAGACAGTGGTCACAGGCTCCCTAGATGACCTGGTGAAGGTCTG
GAAATGGCGTGATGAGAGGCTGGACCTACAGTGGAGTCTGGAGGGACATCAGCTGGGAGTGGTGTCTGTG
GACATCAGCCACACCCTGCCATTGCTGCATCCAGCTCTCTTGATGCTCATATTCGTCTTTGGGACTTGG
AAAATGGCAAACAGATAAAGTCCATAGATGCAGGACCTGTGGATGCCTGGACTTTGGCCTTTTCTCCTGA
TTCCCAGTATCTGGCCACAGGAACCTCATGTCGGGAAAGTGAACATTTTTGGTGTGAAAGTGGGAAAAAG
GAATATCTTTGGACACGAGAGGAAAATTCATTCTTAGTATTGCATATAGTCCTGATGGGAAATACCTAG
CCAGTGGAGCCATAGATGGAATCATCAATATTTTTGATATTGCAACTGGAAAACCTCTGCATACCCTGGA
AGGCCATGCCATGCCATTTCGCTCCTTGACCTTTCCCGGACTCCCAGCTCCTTGCTCACTGCTTCAGAT
GATGGCTACATCAAGATCTATGATGTACAACATGCCAATTTGGCTGGCACGCTGAGCGGCCATGCCTCCT
GGGTGCTGAACGTTGCATTCTGCTCCTGATGACTCACTTTGTTCCAGTTCGTCTGACAAAAGTGTA
AGTTTGGGATGTTGGAACGAGGACTTGTGTTACACCTTCTTTGATCACCAGGATCAGGTCTGGGGAGTA
AAATACAATGAAATGGTTCAAAAATTGTCTGTTGGAGATGACCAGGAAATTCACATCTATGATTGTC
CAATT

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC203963 protein sequence
Red=Cloning site Green=Tags(s)

MTNQYGILFKQEQAHDDAIWSVAWGNTNKKENSETVVTGSLDDLKVKVWWRDERLDLQWSLEGHQLGVVSV
 DISHTLPIAASSSLDAHIRLWDLENGKQIKSIDAGPVDAWTLAFSPDSQYLATGTHVGKVNIFGVESGKK
 EYSLDRGKFILSIAYSPDGKYLASGAIDGIINIFDIATGKLLHTLEGHAMPIRSLTFSPDSQLLVTASD
 DGYIKIYDVQHANLAGTLSGHASWYLNVAFCPDDTHFVSSSSDKSVKVDVGTRTC VHTFFDHQDQVWGV
 KYNGNGSKIVSVGDDQEIHIYDCPI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6422_b10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_025234

ORF Size: 915 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_025234.3](#)

RefSeq Size: 1205 bp

RefSeq ORF: 918 bp

Locus ID: 80349

UniProt ID: [Q9GZS3](#)

Cytogenetics: 15q25.1


Domains: WD40

Protein Pathways: RNA degradation


MW: 33.6 kDa

Gene Summary: WDR61 is a subunit of the human PAF and SKI complexes, which function in transcriptional regulation and are involved in events downstream of RNA synthesis, such as RNA surveillance (Zhu et al., 2005 [PubMed 16024656]).[supplied by OMIM, Mar 2008]


Product images:




Circular map for RC203963



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY WDR61 (Cat# RC203963, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-WDR61 (Cat# [TA807590])(1:2000). Positive lysates [LY410745] (100ug) and [LC410745] (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY410745]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC203963 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified WDR61 protein (Cat# [TP303963]). The protein was produced from HEK293T cells transfected with WDR61 cDNA clone (Cat# RC203963) using MegaTran 2.0 (Cat# [TT210002]).