

Product datasheet for **RC203963L4V**

WDR61 (NM_025234) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WDR61 (NM_025234) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WDR61
Synonyms:	REC14; SKI8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025234
ORF Size:	915 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203963).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025234.1
RefSeq Size:	1205 bp
RefSeq ORF:	918 bp
Locus ID:	80349
UniProt ID:	Q9GZS3
Cytogenetics:	15q25.1
Domains:	WD40
Protein Pathways:	RNA degradation



[View online »](#)

MW: 33.6 kDa

Gene Summary: WDR61 is a subunit of the human PAF and SKI complexes, which function in transcriptional regulation and are involved in events downstream of RNA synthesis, such as RNA surveillance (Zhu et al., 2005 [PubMed 16024656]).[supplied by OMIM, Mar 2008]