

Product datasheet for **RC203964A1V**

Human PLEKHF2 (NM_024613) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PLEKHF2 (NM_024613) AAV Particle
Tag:	Myc-DDK
Symbol:	PLEKHF2
Synonyms:	EAPF; PHAFIN2; ZFYVE18
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203964 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTGGATCGCTTGCCAAACAGTGAAGCAAATACTAGACGTATAAGTATAGTGGAAACTGTTTTGGAG
CAGCTGGTCAACCTTTAACTATACCTGGACGAGTTCTTATTGGAGAAGGAGATTGACTAAGTTGTGCAG
GAAAAAGCCCAAAGCAAGGCAGTTTTTCTTGTTAATGATATTCTTGATATGGCAATATTGTCATCCAG
AAGAAAAATATAACAAACAACATATTATCCCTGGAAAATGTCAGTATTGATCCATCAAAGATGAGG
GAGACTTAAGGAATGGATGGCTAATCAAGACCACTAAATCTTTGCAGTTTATGCTGCCACTGCTAC
GGAGAAATCAGAATGGATGAATCATATAAATAAATGTGTACTGATTTACTCTCCAAAAGTGGGAAGACA
CCCAGTAATGAACATGCTGCTGTCTGGTTCCTGACTCTGAGGCAACTGTATGTATGCGTTGTCAGAAAG
CAAAATTCACACCTGTTAATCGTCGCCACCATGCGCAAAATGTGGTTTTGTTGTCTGTGGGCCCTGCTC
TGAAAAGAGATTTCTTCTCCAGCCAGTCCTCTAAGCCTGTGCGGATTTGTGACTTCTGCTATGACCTG
CTTTCTGCTGGGGACATGGCCACATGCCAGCCTGCTAGATCAGACTCTTACAGCCAGTCATTGAAGTCTC
CTTTAAATGATATGCTGATGATGATGACGATGATGATAGCAGTGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203964 protein sequence Red=Cloning site Green=Tags(s) MYDRLANSEANTRRISIVENCFAAGQPLTIPGRVLIGEGVLTKLCRKKPKARQFFLFNDILVYGNIVIQ KKKYNKQHIIPLENTIDSIDKDEGDLRNGWLITPTKSFAYYAATATEKSEWMNHINKCVTDLLSKSGKT PSNEHAAVWVPDSEATVCMRCQKAKFTPVNRRHHCRCGFVVCGPCSEKRFLFPSQSSKPVRICDFCYDL LSAGDMATCQPARSDSYSQSLKSPLNDMSDDDDDDSSD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024613
ORF Size:	747 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024613.2</u>
RefSeq Size:	2977 bp
RefSeq ORF:	750 bp
Locus ID:	79666
UniProt ID:	<u>Q9H8W4</u>
Cytogenetics:	8q22.1
MW:	27.8 kDa