

## Product datasheet for **RC203965L3V**

### **COPS7A (NM\_016319) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	COPS7A (NM_016319) Human Tagged ORF Clone Lentiviral Particle
Symbol:	COPS7A
Synonyms:	CSN7; CSN7A; SGN7a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016319
ORF Size:	825 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203965).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016319.1</a>
RefSeq Size:	1868 bp
RefSeq ORF:	828 bp
Locus ID:	50813
UniProt ID:	<a href="#">Q9UBW8</a>
Cytogenetics:	12p13.31
Domains:	PCI
MW:	30.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described. [provided by RefSeq, Mar 2014]