

Product datasheet for **RC203967A1V**

Human CHMP6 (NM_024591) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CHMP6 (NM_024591) AAV Particle
Tag: Myc-DDK
Symbol: CHMP6
Synonyms: VPS20
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203967 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTAACTGTTTCGGCCGAAGAAGCAGAGCCGCTCACGGAGCAGGACAAGGCCATCCTGCAACTGA
AGCAGCAGCGGGACAAGCTGAGGCAGTACCAGAAGAGGATCGCCAGCAGCTGGAGCGCAGCGGCCCT
GGCCCGGCAGCTGCTGCGGGACGGCAGGAAGGAACGGGCCAAGCTGCTGCTCAAGAAGAAGCGATACCAG
GAGCAGCTCCTGGACAGGACGGAGAACCAGATCAGCAGCCTGGAGGCCATGGTTCAGAGTATTGAGTTCA
CCCAGATCGAAATGAAAGTGATGGAGGGGCTGCAGTTTGAAATGAGTGTCTGAACAAGATGCACAGGT
GATGTCCATTGAAGAGGTGGAGAGGATCCTGGACGAGACGCAGGAGGCCGTGGAGTACCAGCGGCAATA
GACGAGCTCCTGGCAGGAAGCTTCACTCAGGAGGATGAAGACGCCATCCTGGAGGAGCTGAGCGCAATCA
CTCAGGAACAAATAGAGCTGCCAGAGGTTCCCTCCGAGCCCTTCCTGAGAAGATCCCAGAAAACGTCCC
TGTCAGGCCAGGCCAGGCAGGCGGAGCTGGTGGCAGCTTCG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203967 protein sequence
Red=Cloning site Green=Tags(s)

MGNLFGRKKQSRVTEQDKAILQLKQQRDKLRQYQKRIAQQLERERALARQLLRDGRKERAKLLKKKRYQ
EQLLDRTENQISSLEAMVQSIEFTQIEMKVMGLQFGNECLNMHVMSIEEVERILDETQEAIVEYQRQI
DELLAGSFTQEDEDAILEELSAITQEIQIELPEVPSEPLPEKIPENVVKARPRQAEVAAS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_024591
ORF Size:	603 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024591.3</u>
RefSeq Size:	1678 bp
RefSeq ORF:	606 bp
Locus ID:	79643
UniProt ID:	<u>Q96FZ7</u>
Cytogenetics:	17q25.3
MW:	23.5 kDa