

Product datasheet for **RC203968L4V**

BEN Domain Containing Protein 5 (BEND5) (NM_024603) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BEN Domain Containing Protein 5 (BEND5) (NM_024603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BEND5
Synonyms:	C1orf165
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024603
ORF Size:	1263 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203968).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024603.3
RefSeq Size:	1781 bp
RefSeq ORF:	1266 bp
Locus ID:	79656
UniProt ID:	Q7L4P6
Cytogenetics:	1p33
MW:	48.2 kDa



[View online »](#)

Gene Summary: Acts as a transcriptional repressor (PubMed:23468431).[UniProtKB/Swiss-Prot Function]