

## Product datasheet for **RC203971L4V**

### **APOBEC3C (NM\_014508) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	APOBEC3C (NM_014508) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APOBEC3C
Synonyms:	A3C; APOBEC1L; ARDC2; ARDC4; ARP5; bK150C2.3; PBI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014508
ORF Size:	570 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC203971).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014508.2</a>
RefSeq Size:	1127 bp
RefSeq ORF:	573 bp
Locus ID:	27350
UniProt ID:	<a href="#">Q9NRW3</a>
Cytogenetics:	22q13.1
MW:	22.8 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control. [provided by RefSeq, Jul 2008]