

Product datasheet for **RC203990A1V**

Human RCN1 (NM_002901) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RCN1 (NM_002901) AAV Particle
Tag:	Myc-DDK
Symbol:	RCN1
Synonyms:	HEL-S-84; PIG20; RCAL; RCN
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC203990 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGCGCGGTGGCGCGGCCCGCCCTGGGGTTAGCCCTGGGGCTGCTGCTGGCGCTGGTCTGGCGC
CGCGGGTTCTGCGGGCCAAGCCACGGTGCACAAAGAGCGCGTGGTGCGGCCCGACTCGGAGCTGGGCGA
GCGGCCCCCTGAGGACAACAGAGCTTCCAGTACGACCACGAGGCCTTCCTGGCAAGGAGGACTCCAAG
ACCTTCGACAGCTCACCCCGGACGAGAGCAAGGAGAGGCTAGGAAGATTGTTGATCGAATCGACAATG
ATGGGGATGGCTTTGTCACTACTGAGGAGCTGAAAACCTGGATCAAACGGGTGCAGAAAAGATACATCTT
TGATAATGTGCGCAAAGTCTGGAAGGATTATGATAGGGACAAGGATGATAAAATTCCTGGGAAGAATAC
AAACAAGCCACCTATGGTTACTACCTAGGAAACCCGACAGATTTTCATGATTCTTCAGATCATCACCT
TTAAAAAGATGCTGCCACGTGATGAGAGAAGATTCAAAGCTGCAGACCTCAATGGTGACCTGACAGCTAC
TCGGGAGGAGTTCACTGCCTTTCTGCATCCTGAAGAGTTTGAACATATGAAGGAAATTGTGGTTTTGAA
ACCCTGGAGGACATCGACAAGAACGGGGATGGGTTTGTGGATCAGGATGAGTATATTGCGGATATGTTTT
CCCATGAGGAGAATGGCCCTGAGCCAGACTGGGTTTTATCAGAACGGGAGCAGTTTAACGAATTCGGGA
TCTGAACAAGGACGGGAAGTTAGACAAAGATGAGATTCGCCACTGGATCCTCCCTCAAGATTATGATCAT
GCACAGGCTGAGGCCAGGCATCTGGTATATGAATCAGACAAAAACAAGGATGAGAAGCTAACTAAAGAGG
AAATATTGGAGAACTGGAACATGTTTGTGCGGAAGCCAAGCTACCAATTACGGGGAAGATCTCACAAAAA
TCATGATGAGCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC203990 protein sequence Red=Cloning site Green=Tags(s) MARGGRGRRLLGLALGLLLALVLAPRVLRAKPTVRKERVVRPDSELGERPPEDNQSFQYDHEAFLGKEDSK TFDQLTPDESKERLGKIVDRIDNDGDFVTTEELKTWIKRVQKRYIFDNVAKVWKDYDRDKDDKISWEEY KQATYGYLLGNPAEFHDSSDHHTFKKMLPRDERRFKAADLNGDLTATREEFTAFLEHPEEFHEMKEIVVLE TLEDIDKNGDGFVDQDEYIADMFSHEENGPEPDWVLSEREQFNEFRDLNKDGKLDKDEIRHWILPQDYDH AQAEARHLVYESDKNKDEKLTKEEILENWNMFVGSQATNYGEDLTKNHDEL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002901
ORF Size:	993 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002901.1</u>
RefSeq Size:	2529 bp
RefSeq ORF:	996 bp
Locus ID:	5954
UniProt ID:	<u>Q15293</u>
Cytogenetics:	11p13
MW:	38.9 kDa