

Product datasheet for **RC203995A1V**

Human Midkine (MDK) (NM_001012333) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Midkine (MDK) (NM_001012333) AAV Particle
Tag: Myc-DDK
Symbol: Midkine
Synonyms: ARAP; MK; NEGF2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC203995 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCAGCACCGAGGCTTCCTCCTCACCCCTCTCGCCCTGCTGGCGCTCACCTCCGCGGTGCGCAAAA
 AGAAAGATAAGGTGAAGAAGGGCGGCCCGGGGAGCGAGTGCCTGAGTGGGCCTGGGGCCCTGCACCCC
 CAGCAGCAAGGATTGCGGCGTGGGTTTCCGCGAGGGCACCTGCGGGGCCAGACCCAGCGCATCCGGTGC
 AGGGTGCCTGCAACTGGAAGAAGGAGTTTGGAGCCGACTGCAAGTACAAGTTTGAGAACTGGGGTGCCT
 GTGATGGGGGCACAGGCACCAAAGTCCGCCAAGGCACCCTGAAGAAGGGCGCTACAATGCTCAGTGCCA
 GGAGACCATCCGCGTCACCAAGCCCTGCACCCCAAGACCAAAGGCAAGGCCAAAGCCAAGAAAGGGAAG
 GGAAAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC203995 protein sequence
 Red=Cloning site Green=Tags(s)

MQHRGFLLLTLALLALTSAVAKKKDKVKGGPGSECAEWAWGPCTPSSKDCGVGFREGTCGAQTQRIRC
 RVPNCNWKKEFGADCKYKFENWGACDGGTGKVRQGTLLKARYNAQCQETIRVTKPCTPKTKAKAKKGGK
 GKD

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_001012333
ORF Size:	429 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001012333.1</u>
RefSeq Size:	969 bp
RefSeq ORF:	432 bp
Locus ID:	4192
UniProt ID:	<u>P21741</u>
Cytogenetics:	11p11.2
MW:	15.6 kDa