

Product datasheet for **RC204002A1V**

Human TST (NM_003312) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TST (NM_003312) AAV Particle
Tag:	Myc-DDK
Symbol:	TST
Synonyms:	RDS
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204002 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTTCATCAGGTGCTCTACCGGGCGCTGGTCTCCACCAAGTGGCTGGCGGAGTCCATCAGGACTGGCA
AGCTGGGGCCCGCCTGCGGGTCTGGACGCGTCCTGGTACTCACCAGGCACCCGAGAGGCCGCAAGGA
GTACCTCGAGCGCCACGTACCGGGCGCCTCTTTCTTTGACATAGAAGAGTGCCGGGACACGGCGTCGCCC
TACGAGATGATGCTGCCAGCGAGGCTGGCTTCGCCGAGTATGTGGGCCGCTGGGCATCAGCAACCACA
CGCACGTGGTGGTGTATGATGGTGAACACCTGGGCAGTCTCTATGCTCCCGGGTCTGGTGGATGTTCCG
TGTGTTTGGCCACCGCACCGTATCAGTGCTCAATGGTGGCTTCCGGAAGTGGCTGAAGGAGGGCCACCCG
GTGACATCCGAGCCCTCACGCCAGAACCGGCCGTCTTCAAAGCCACACTGGACCGCTCCCTGCTCAAGA
CCTACGAGCAGGTGCTGGAGAACCTTGAATCTAAGAGGTTCCAGCTGGTGGATTCAAGGTCTCAAGGGCG
GTTCTGGGCACCGAGCCGAGCCGGATGCAGTAGGACTGGACTCGGGCCATATCCGTGGTGCCGTCAAC
ATGCCTTTCATGGACTTCCTGACTGAGGATGGCTTCGAGAAGGGCCAGAAAGAGTCCGTGCTCTGTTCC
AGACCAAGAAGGTGGATCTCTCGCAGCCTCTATTGCCACGTGCCGAAGGGAGTCAACCGCTGCCACGT
GGCCTTGGCTGCCTACCTCTGCGCAAGCCTGATGTGGCCGTGTACGATGGCTCCTGGTCCGAGTGGTTT
CGCCGGGCCCCCAGAGAGCCGTGTGTCCAGGAAAGTCTGAGAAGGCC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204002 protein sequence Red=Cloning site Green=Tags(s) MVHQVLRYALVSTKWLAEISRTGKLGPLRVLDASWYSPGTREARKEYLERHVPASFFDIEECRDTASP YEMMLPSEAGFAEYVGRLGISNHTHVVVYDGEHLGSFYAPRVWWMFRVFGHRTVSVLNGGFRNWLKEGHP VTSEPSRPEPAVFKATLDRSLKTYEQVLENLESKRFQLVDSRSQGRFLGTEPEPDAVGLDSGHIRGAVN MPFMDFLTEDGFEGPEELRALFQTKKVDLSQPLIATCRKGVTAACHVALAAYLCGKPDVAVYDGSWSEWF RRAPPESRVSQGKSEKA SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003312
ORF Size:	891 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003312.4</u>
RefSeq Size:	1143 bp
RefSeq ORF:	894 bp
Locus ID:	7263
UniProt ID:	<u>Q16762</u>
Cytogenetics:	22q12.3
MW:	33.4 kDa