

Product datasheet for **RC204006A1V**

Human RPL11 (NM_000975) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RPL11 (NM_000975) AAV Particle
Tag: Myc-DDK
Symbol: RPL11
Synonyms: DBA7; GIG34; L11; uL5
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204006 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGATCAAGGTGAAAAGGAGAACCCCATGCGGGAACCTTCGCATCCGCAAACCTGTCTCAACATCT
 GTGTTGGGGAGAGTGGAGACAGACTGACGCGAGCAGCCAAGGTGTTGGAGCAGCTCACAGGGCAGACCCC
 TGTGTTTTCCAAAGCTAGATACACTGTCAGATCCTTTGGCATCCGGAGAAATGAAAAGATTGCTGTCCAC
 TGCACAGTTCGAGGGGCCAAGGCAGAGAAATCTTGAGAAGGGTCTAAAGGTGCGGGAGTATGAGTTAA
 GAAAAACAACCTTCTCAGATACTGGAACCTTTGGTTTTGGGATCCAGGAACACATCGATCTGGGTATCAA
 ATATGACCCAAGCATTGGTATCTACGGCCTGGACTTCTATGTGGTGCTGGGTAGGCCAGGTTTCAGCATC
 GCAGACAAGAAGCGCAGGACAGGCTGCATTGGGGCCAAACACAGAATCAGCAAAGAGGAGGCCATGCGCT
 GGTTCAGCAGAAGTATGATGGGATCATCCTTCTGGCAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204006 protein sequence
 Red=Cloning site Green=Tags(s)

MADQGEKENPMRELRIKRLCLNICVGESGDRL TRAAKVLEQL TGQTPVFSKARYTVRSFGIRRNEKIAVH
 CTVRGAKAEIELEKLVREYELRKNFSDTGNFGFGIQEHIDLGIKYDPSIGIYGLDFYVVLGRPFGSI
 ADKKRRITGCIGAKHRISKEEAMRWFQQKYDGIILPGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_000975
ORF Size:	531 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_000975.2
RefSeq Size:	644 bp
RefSeq ORF:	537 bp
Locus ID:	6135
UniProt ID:	P62913
Cytogenetics:	1p36.11
MW:	20.1 kDa