

## Product datasheet for **RC204007L4V**

### Tropomyosin 2 (TPM2) (NM\_213674) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tropomyosin 2 (TPM2) (NM_213674) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Tropomyosin 2
Synonyms:	AMCD1; DA1; DA2B; DA2B4; HEL-S-273; NEM4; TMSB
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_213674
ORF Size:	852 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204007).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_213674.1</a> , <a href="#">NP_998839.1</a>
RefSeq Size:	1182 bp
RefSeq ORF:	855 bp
Locus ID:	7169
UniProt ID:	<a href="#">P07951</a>
Cytogenetics:	9p13.3
Protein Pathways:	Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM)
MW:	33 kDa



[View online »](#)

**Gene Summary:**

This gene encodes beta-tropomyosin, a member of the actin filament binding protein family, and mainly expressed in slow, type 1 muscle fibers. Mutations in this gene can alter the expression of other sarcomeric tropomyosin proteins, and cause cap disease, nemaline myopathy and distal arthrogyriposis syndromes. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Mar 2009]