

## Product datasheet for **RC204009L3V**

### **HARS (HARS1) (NM\_002109) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HARS (HARS1) (NM_002109) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HARS1
Synonyms:	CMT2W; HARS; HRS; USH3B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002109
ORF Size:	1527 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204009).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002109.3</a>
RefSeq Size:	1981 bp
RefSeq ORF:	1530 bp
Locus ID:	3035
UniProt ID:	<a href="#">P12081</a>
Cytogenetics:	5q31.3
Domains:	WHEP-TRS, tRNA-synt_2b, HGTP_anticodon
Protein Pathways:	Aminoacyl-tRNA biosynthesis



[View online »](#)

MW: 57.2 kDa

**Gene Summary:** Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate amino acids. The protein encoded by this gene is a cytoplasmic enzyme which belongs to the class II family of aminoacyl-tRNA synthetases. The enzyme is responsible for the synthesis of histidyl-transfer RNA, which is essential for the incorporation of histidine into proteins. The gene is located in a head-to-head orientation with HARSL on chromosome five, where the homologous genes share a bidirectional promoter. The gene product is a frequent target of autoantibodies in the human autoimmune disease polymyositis/dermatomyositis. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]