

Product datasheet for **RC204011L3V**

SF20 (MYDGF) (NM_019107) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SF20 (MYDGF) (NM_019107) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SF20
Synonyms:	C19orf10; EUROIMAGE1875335; IL25; IL27; IL27w; R33729_1; SF20
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019107
ORF Size:	519 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204011).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019107.3
RefSeq Size:	1067 bp
RefSeq ORF:	522 bp
Locus ID:	56005
UniProt ID:	Q969H8
Cytogenetics:	19p13.3
Protein Families:	Secreted Protein
MW:	18.8 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene was previously thought to support proliferation of lymphoid cells and was considered an interleukin. However, this activity has not been reproducible and the function of this protein is currently unknown. [provided by RefSeq, Jul 2008]