

Product datasheet for **RC204018L4V**

MAGEF1 (NM_022149) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MAGEF1 (NM_022149) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAGEF1
Synonyms:	MAGE-F1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022149
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204018).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022149.4
RefSeq Size:	1700 bp
RefSeq ORF:	924 bp
Locus ID:	64110
UniProt ID:	Q9HAY2
Cytogenetics:	3q27.1
Domains:	MAGE
MW:	35.3 kDa


[View online »](#)

Gene Summary:

This intronless gene encodes a member of the MAGE superfamily. It is ubiquitously expressed in normal tissues and in tumor cells. This gene includes a microsatellite repeat in the coding region. [provided by RefSeq, Jul 2008]