

## Product datasheet for **RC204023L2V**

### TRAX (TSNAX) (NM\_005999) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRAX (TSNAX) (NM_005999) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRAX
Synonyms:	C3PO; TRAX
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005999
ORF Size:	870 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204023).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005999.2</a>
RefSeq Size:	2667 bp
RefSeq ORF:	873 bp
Locus ID:	7257
UniProt ID:	<a href="#">Q99598</a>
Cytogenetics:	1q42.2
Domains:	Translin
MW:	33.1 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein which specifically interacts with translin, a DNA-binding protein that binds consensus sequences at breakpoint junctions of chromosomal translocations. The encoded protein contains bipartite nuclear targeting sequences that may provide nuclear transport for translin, which lacks any nuclear targeting motifs. [provided by RefSeq, Jul 2008]