

Product datasheet for **RC204026A1V**

Human MED8 (NM_201542) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human MED8 (NM_201542) AAV Particle
Tag: Myc-DDK
Symbol: MED8
Synonyms: ARC32
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204026 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCGGCAGACTGAAGGACGGGTGCCTGTTTTAGCCATGAGGTAGTCCCTGACCATCTGAGAACCAAGC
 CTGACCTGAAGTGGAAGAAGCAGGAGAAGCACTGACGACAGATGCTGCCCGCATTGGTGCAGATGCAGC
 CCAGAAGCAGATCCAGAGCTTGAATAAAATGTGTTCAAACCTTCTGGAGAAAAATCAGCAAAGAGGAGCGA
 GAATCAGAGAGTGGAGGTCTCCGGCCGAACAAGCAGACCTTTAACCTACAGACCTAATGCCTTGGTGG
 CAGCTGTTGCCTTTGGAAAGGACTATCTAATTGGAGACCTTCAGGCAGCAGTGGTCTGGCCAGGCAGG
 CAGCCAGGAGCTGGGACGATCCTTGCAGGAACCTCAGGATTACAGCAGGTGCAGATGGCAGGAGCTCCA
 AGCCAGCAGCAGCCAATGCTCAGTGGGGTACAAATGGCTCAGGCAGGTCAACCAGGAAAATGCCAAGTG
 GAATAAAACCAACATCAAGTCGGCTTCATGCATCCCTACCAGCGG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204026 protein sequence
 Red=Cloning site Green=Tags(s)

MRQTEGRVPVFSHEVVPDHLRTKPDPEVEEQEKQLTTDAARIGADAAQKQIQSLNKMCSNLLEKISKEER
 ESESGGLRPNKQTFNPTDTNALVAAVAFGKLSNWRPSGSSGPGQAGQPGAGTILAGTSGLQQVQMAGAP
 SQQPMLSGVQMAQAGQPGKMPSGIKTNIKSASMPHYQR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_201542
ORF Size:	537 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_201542.3</u> , <u>NP_963836.1</u>
RefSeq Size:	1989 bp
RefSeq ORF:	807 bp
Locus ID:	112950
UniProt ID:	<u>Q96G25</u>
Cytogenetics:	1p34.2
MW:	19 kDa