

## Product datasheet for **RC204026L3V**

### MED8 (NM\_201542) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MED8 (NM_201542) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MED8
Synonyms:	ARC32
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_201542
ORF Size:	537 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204026).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_201542.3</a> , <a href="#">NP_963836.1</a>
RefSeq Size:	1989 bp
RefSeq ORF:	807 bp
Locus ID:	112950
UniProt ID:	<a href="#">Q96G25</a>
Cytogenetics:	1p34.2
Protein Families:	Druggable Genome, Transcription Factors
MW:	19 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein component of the mediator complex, which aids in transcriptional activation through interaction with RNA polymerase II and gene-specific transcription factors. The encoded protein may also function in ubiquitin ligation and protein degradation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2013]