

## Product datasheet for **RC204029L4V**

### **SAMM50 (NM\_015380) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

<b>Product Type:</b>	Lentiviral Particles
<b>Product Name:</b>	SAMM50 (NM_015380) Human Tagged ORF Clone Lentiviral Particle
<b>Symbol:</b>	SAMM50
<b>Synonyms:</b>	CGI-51; OMP85; SAM50; TOB55; TRG-3; YNL026W
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Vector:</b>	pLenti-C-mGFP-P2A-Puro (PS100093)
<b>Tag:</b>	mGFP
<b>ACCN:</b>	NM_015380
<b>ORF Size:</b>	1407 bp
<b>ORF Nucleotide Sequence:</b>	The ORF insert of this clone is exactly the same as(RC204029).
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>RefSeq:</b>	<a href="#">NM_015380.3</a> , <a href="#">NP_056195.2</a>
<b>RefSeq Size:</b>	1773 bp
<b>RefSeq ORF:</b>	1410 bp
<b>Locus ID:</b>	25813
<b>UniProt ID:</b>	<a href="#">Q9Y512</a>
<b>Cytogenetics:</b>	22q13.31
<b>Domains:</b>	Bac_surface_Ag
<b>MW:</b>	52 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a component of the Sorting and Assembly Machinery (SAM) of the mitochondrial outer membrane. The Sam complex functions in the assembly of beta-barrel proteins into the outer mitochondrial membrane.[provided by RefSeq, Jun 2011]