

Product datasheet for **RC204030**

C9orf16 (NM_024112) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C9orf16 (NM_024112) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C9orf16
Synonyms: EST00098
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC204030 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCGGGCCCAACGGAGACCTGGGGATGCCGGTGGAGCGGGAGCGGAAGCGGAGGAGGACGGCTTCG
GGGAAGCAGAATACGCTGCCATCAACTCCATGCTGGACCAGATCAACTCCTGTCTGGACCACCTGGAGGA
GAAGAATGACCACCTCCACGCCCGCTCCAGGAGCTGCTGGAGTCCAACCGGCAGACACGCCCTGGAGTTC
CAGCAGCAGCTCGGGGAGCCCCAGTGATGCCAGCCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204030 protein sequence
Red=Cloning site Green=Tags(s)
MSGPNGDLGMPVEAGAEGEEDGFGEAEYAAINSMLDQINSCLDHLEEKNDHLHARLQELLESNRQTRLEF
QQQLGEAPSDASP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6422_d05.zip

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_024112

ORF Size: 249 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024112.4](#)


RefSeq Size: 723 bp

RefSeq ORF: 252 bp


Locus ID: 79095

UniProt ID: Q9BUW7
 Cytogenetics: 9q34.11
 Domains: UPF0184
 MW: 9.1 kDa


Product images:



Circular map for RC204030



Western blot validation of overexpression lysate (Cat# [LY411346]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC204030 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified C9orf16 protein (Cat# [TP304030]). The protein was produced from HEK293T cells transfected with C9orf16 cDNA clone (Cat# RC204030) using MegaTran 2.0 (Cat# [TT210002]).