

Product datasheet for **RC204030A1V**

Human BBLN (NM_024112) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human BBLN (NM_024112) AAV Particle
Tag: Myc-DDK
Symbol: BBLN
Synonyms: EST00098
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204030 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCGGGCCCAACGGAGACCTGGGGATGCCGGTGGAGGCGGGAGCGGAAGGCGAGGAGGACGGCTTCG
GGGAAGCAGAAATACGCTGCCATCAACTCCATGCTGGACCAGATCAACTCCTGTCTGGACCACCTGGAGGA
GAAGAATGACCACCTCCACGCCCGCCTCCAGGAGCTGCTGGAGTCCAACCGGCAGACACGCCTGGAGTTC
CAGCAGCAGCTCGGGGAGGCCCCAGTGATGCCAGCCCC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204030 protein sequence
Red=Cloning site Green=Tags(s)

MSGPNGDLGMPVEAGAEGEEDGFGEAEYAAINSMLDQINSCLDHLEEKNDHLHARLQELLESNRQTRLEF
QQQLGEAPSDASP

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_024112
ORF Size: 249 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024112.2</u>
RefSeq Size:	723 bp
RefSeq ORF:	252 bp
Locus ID:	79095
UniProt ID:	<u>Q9BUW7</u>
Cytogenetics:	9q34.11
MW:	9.1 kDa