

Product datasheet for **RC204060L3V**

IGFBP6 (NM_002178) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IGFBP6 (NM_002178) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IGFBP6
Synonyms:	IBP6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002178
ORF Size:	720 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204060).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002178.2
RefSeq Size:	980 bp
RefSeq ORF:	723 bp
Locus ID:	3489
UniProt ID:	P24592
Cytogenetics:	12q13.13
Protein Families:	Secreted Protein
MW:	25.3 kDa



[View online »](#)

Gene Summary:

IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.[UniProtKB/Swiss-Prot Function]