

Product datasheet for **RC204065L4V**

TTC17 (NM_018259) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TTC17 (NM_018259) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TTC17
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018259
ORF Size:	3423 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204065).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018259.3
RefSeq Size:	4525 bp
RefSeq ORF:	3426 bp
Locus ID:	55761
UniProt ID:	Q96AE7
Cytogenetics:	11p12-p11.2
Domains:	TPR
Protein Families:	Protease
MW:	129.6 kDa



[View online »](#)

Gene Summary: Plays a role in primary ciliogenesis by modulating actin polymerization.[UniProtKB/Swiss-Prot Function]