

## Product datasheet for **RC204067L4V**

### EXT2 (NM\_207122) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EXT2 (NM_207122) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EXT2
Synonyms:	SOTV; SSMS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207122
ORF Size:	2154 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204067).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207122.1</a>
RefSeq Size:	3651 bp
RefSeq ORF:	2157 bp
Locus ID:	2132
UniProt ID:	<a href="#">Q93063</a>
Cytogenetics:	11p11.2
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Heparan sulfate biosynthesis, Metabolic pathways



[View online »](#)

**MW:** 82.1 kDa

**Gene Summary:** This gene encodes one of two glycosyltransferases involved in the chain elongation step of heparan sulfate biosynthesis. Mutations in this gene cause the type II form of multiple exostoses. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Jul 2008]