

Product datasheet for **RC204068L2V**

p66 alpha (GATAD2A) (NM_017660) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | p66 alpha (GATAD2A) (NM_017660) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | p66 alpha |
| Synonyms: | p66alpha |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_017660 |
| ORF Size: | 1470 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC204068). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_017660.3 , NP_060130.2 |
| RefSeq Size: | 5689 bp |
| RefSeq ORF: | 1902 bp |
| Locus ID: | 54815 |
| UniProt ID: | Q86YP4 |
| Cytogenetics: | 19p13.11 |
| Protein Families: | Transcription Factors |
| MW: | 52.5 kDa |



[View online »](#)

Gene Summary:

Transcriptional repressor. Enhances MBD2-mediated repression. Efficient repression requires the presence of GATAD2B.[UniProtKB/Swiss-Prot Function]