

Product datasheet for **RC204068L3V**

p66 alpha (GATAD2A) (NM_017660) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	p66 alpha (GATAD2A) (NM_017660) Human Tagged ORF Clone Lentiviral Particle
Symbol:	p66 alpha
Synonyms:	p66alpha
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017660
ORF Size:	1470 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204068).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017660.3 , NP_060130.2
RefSeq Size:	5689 bp
RefSeq ORF:	1902 bp
Locus ID:	54815
UniProt ID:	Q86YP4
Cytogenetics:	19p13.11
Protein Families:	Transcription Factors
MW:	52.5 kDa



[View online »](#)

Gene Summary:

Transcriptional repressor. Enhances MBD2-mediated repression. Efficient repression requires the presence of GATAD2B.[UniProtKB/Swiss-Prot Function]