

Product datasheet for **RC204070L4V**

SFT2D1 (NM_145169) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SFT2D1 (NM_145169) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SFT2D1
Synonyms:	C6orf83; pRGR1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145169
ORF Size:	477 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204070).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145169.1
RefSeq Size:	717 bp
RefSeq ORF:	480 bp
Locus ID:	113402
UniProt ID:	Q8WV19
Cytogenetics:	6q27
Protein Families:	Transmembrane
MW:	17.8 kDa



[View online »](#)

Gene Summary:

May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.[UniProtKB/Swiss-Prot Function]