

## Product datasheet for **RC204071A1V**

### Human C14orf172 (TRMT61A) (NM\_152307) AAV Particle

#### Product data:

|                           |                                                                   |
|---------------------------|-------------------------------------------------------------------|
| Product Type:             | AAV Particles                                                     |
| Product Name:             | Human C14orf172 (TRMT61A) (NM_152307) AAV Particle                |
| Tag:                      | Myc-DDK                                                           |
| Symbol:                   | C14orf172                                                         |
| Synonyms:                 | C14orf172; GCD14; Gcd14p; hTRM61; TRM61                           |
| Mammalian Cell Selection: | None                                                              |
| Vector:                   | pAAV-AC-Myc-DDK (PS100089)                                        |
| ORF Nucleotide Sequence:  | >RC204071 ORF sequence<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGCTTCGTGGCATAACGAGGAGCTGATCAAGGAGGGTGACACGGCCATCCTGTCACTGGGCCATGGTG  
CAATGGTGGCAGTGCGTGTGCAGCGTGGGGCACAGACCCAGACCCGGCATGGTGTCTGCGGCACTCAGT  
TGACCTTATCGGCCGCCCTTCGGCTCCAAGGTGACGTGCGGCCGAGGTGGTGGGTGTATGTGCTGCAC  
CCCACGCCCGAGCTCTGGACGCTGAACCTGCCGACCGCACGCAGATCCTCTACTCCACAGACATCGCCC  
TCATCACCATGATGTTGGAGCTTCGGCCCGGCTCTGTGGTCTGTGAGTCTGGCACCAGGAGTGGCTCTGT  
GTCCACGCCATCATCCGACCAATTGCACCCACGGGTACCTGCACACGGTGGAGTTCCACAGCAGCGG  
GCAGAGAAGGCCCGGGAGGAGTCCAGGAGCACCGTGTGGGCCGTGGGTGACTGTGCGCACCCAGGACG  
TGTGCCGAGTGGCTTTGGCGTGAGCCACGTGGCCGACGCCGTCTTCTGGACATCCCATCACCTGGGA  
GGCCGTGGGCCACGCTGGGACGCCCTCAAGGTCGAAGGCGGGCGCTTCTGCTCCTTCTACCGTGCATC  
GAGCAGGTGAACGCACATGCCAGGCGTGGCAGCGCGCGCTTCTCAGAGCTGAGCACCTGGAGGTGC  
TGCCACAGGTCTACAACGTGCGCACTGTAGCCTGCCACCGCCGACCTGGGCACAGGCACAGATGGCCC  
TGCCGGCTCCGACACAGCCCTTCCGAGCGGCACGCCATGAAGGAGGCCGTGGGCCACACCGGCTAC  
CTGACCTTCGCCACCAAGACCCAGGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Protein Sequence: | <p>&gt;RC204071 protein sequence</p> <p>Red=Cloning site Green=Tags(s)</p> <p>MSFVAYEELIKEGDTAILSLGHGAMVAVRVQRGAQTQTRHGVLRHVSVDLIGRPFGSKVTCGRGGWVYVLH<br/> PTPELWTLNLPHTQILYSTDIALITMMLELRPGSVVCESGTGSGSVSHAIIRTIAPTGHLHTVEFHQQR<br/> AEKAREEFQEHVRVGRWTVRTQDVCRSFGVSHVADAVFLDIPSPWEAVGHAWDALKVEGGRFCFSFPCI<br/> EQVQRTCQALAARGFSELSTLEVLPQVYNVRTVSLPPDLGTGTDGPAGSDTSPFRSGTPMKEAVGHTGY<br/> LTFATKTPG</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>                            |
| Species:          | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Serotype:         | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| ACCN:             | NM_152307                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| ORF Size:         | 867 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Buffer:           | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Stability:        | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| RefSeq:           | <a href="#">NM_152307.2</a> , <a href="#">NP_689520.2</a>                                                                                                                                                                                                                                                                                                                                                                                                                   |
| RefSeq Size:      | 3257 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| RefSeq ORF:       | 870 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Locus ID:         | 115708                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| UniProt ID:       | <a href="#">Q96FX7</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Cytogenetics:     | 14q32.33                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| MW:               | 31.4 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |