

Product datasheet for **RC204073L1V**

Nogo Receptor (RTN4R) (NM_023004) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Nogo Receptor (RTN4R) (NM_023004) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Nogo Receptor |
| Synonyms: | NGR; NOGOR |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_023004 |
| ORF Size: | 1419 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC204073). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_023004.5 |
| RefSeq Size: | 1924 bp |
| RefSeq ORF: | 1422 bp |
| Locus ID: | 65078 |
| UniProt ID: | Q9BZR6 |
| Cytogenetics: | 22q11.21 |
| Domains: | LRRCT, LRR, LRR_TYP |
| Protein Families: | Druggable Genome |



[View online »](#)

MW: 50.7 kDa

Gene Summary: This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system. [provided by RefSeq, Jul 2008]