

## Product datasheet for **RC204084L2V**

### ERGIC3 (NM\_015966) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ERGIC3 (NM_015966) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ERGIC3
Synonyms:	C2orf47; C20orf47; CGI-54; dj477O4.2; Erv46; NY-BR-84; PRO0989; SDBCAG84
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015966
ORF Size:	1149 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204084).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015966.2</a>
RefSeq Size:	1368 bp
RefSeq ORF:	1152 bp
Locus ID:	51614
UniProt ID:	<a href="#">Q9Y282</a>
Cytogenetics:	20q11.22
Protein Families:	Transmembrane
MW:	43.2 kDa



[View online »](#)

**Gene Summary:** Possible role in transport between endoplasmic reticulum and Golgi.[UniProtKB/Swiss-Prot Function]