

## Product datasheet for **RC204085A1V**

### Human PPCDC (NM\_021823) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human PPCDC (NM\_021823) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** PPCDC  
**Synonyms:** coaC; MDS018; PPC-DC  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC204085 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGAACCAAAGGCCTCCTGTCCAGCTGCTGCACCCTTGATGGAGAGAAAATTCATGTTCTTGTTGGGTG  
TCACGGGGAGTGTGCGAGCCCTGAAGTTGCCTCTTCTGGTGTCAAAGCTTTTGGACATTCCTGGGCTGGA  
AGTAGCAGTGGTCACAACCTGAGAGAGCCAAACATTTCTACAGCCCCAGGACATTCCTGTCAACCTCTAC  
AGCGACGCTGATGAATGGGAGATGTGGAAGAGCCGCTCTGACCCAGTTCTGCACATTGACCTGCGGAGGT  
GGCGACACCTCCTGCTGGTGGCTCCTCTTGATGCCAACACTCTGGGGAAGGTGGCCAGTGGCATCTGTGA  
CAACTTGCTTACCTGCGTCATGCGGGCCTGGGACCGCAGCAAGCCCTGCTCTTCTGCCCCGCCATGAAC  
ACCGCCATGTGGGAGACCCGATCACAGCGCAGCAGGTAGACCAGCTCAAGGCCTTTGGCTATGTCGAGA  
TCCCCTGTGTGGCCAAGAAGCTGGTGTGCGGAGATGAAGGTCTCGGGGCCATGGCTGAAGTGGGACCAT  
CGTGGACAAAGTGAAAGAAGTCCTCTTCAGCACAGTGGCTTCCAGCAGAGT

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC204085 protein sequence  
Red=Cloning site Green=Tags(s)

MEPKASCPAAAPLMERKFHVLVGVTGSVAALKLPLLVSLLDIPGLEAVVTTTERAKHFYSPQDIPVTLY  
SDADEWEMWKS RSDPVLHIDLRRWADLLL VAPLDANTLGKVASGICDNLCTVMRAWDRSKPLFCPAMN  
TAMWEHPITAQQVDQLKAFGYVEIPCVAKKL VCGDEGLGMAEVGTIVDKVKEVL FQHS GFQQS

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human



[View online »](#)

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Serotype:</b>     | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>         | NM_021823                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>ORF Size:</b>     | 612 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <u><a href="#">NM_021823.3</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq Size:</b>  | 2268 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>   | 615 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Locus ID:</b>     | 60490                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>UniProt ID:</b>   | <u><a href="#">Q96CD2</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Cytogenetics:</b> | 15q24.2                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>MW:</b>           | 22.4 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |