

Product datasheet for **RC204113L3V**

RNMTL1 (MRM3) (NM_018146) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RNMTL1 (MRM3) (NM_018146) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNMTL1
Synonyms:	RMTL1; RNMTL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018146
ORF Size:	1260 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204113).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018146.2
RefSeq Size:	1815 bp
RefSeq ORF:	1263 bp
Locus ID:	55178
UniProt ID:	Q9HC36
Cytogenetics:	17p13.3
Domains:	SpoU_methylase
Protein Families:	Stem cell - Pluripotency



[View online »](#)

MW: 47 kDa

Gene Summary: Efficient translation of mitochondrial-derived transcripts requires proper assembly of the large subunit of the mitochondrial ribosome. Central to the biogenesis of this large subunit is the A-loop of mitochondrial 16S rRNA, which is modified by three rRNA methyltransferases located near mtDNA nucleoids. The protein encoded by this gene methylates G(1370) of 16S rRNA, and this modification is necessary for proper ribosomal large subunit assembly. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2015]