

Product datasheet for **RC204117L4V**

RNF31 (NM_017999) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RNF31 (NM_017999) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RNF31
Synonyms:	HOIP; Paul; ZIBRA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017999
ORF Size:	3216 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204117).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017999.4
RefSeq Size:	3571 bp
RefSeq ORF:	3219 bp
Locus ID:	55072
UniProt ID:	Q96EP0
Cytogenetics:	14q12
Domains:	IBR
Protein Families:	Druggable Genome



[View online »](#)

MW: 119.5 kDa

Gene Summary: The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. The encoded protein is the E3 ubiquitin-protein ligase component of the linear ubiquitin chain assembly complex. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2015]