

Product datasheet for RC204131L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

SLC33A1 (NM_004733) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: SLC33A1 (NM 004733) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: SLC33A1

Synonyms: ACATN; AT-1; AT1; CCHLND; SPG42

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC204131).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_004733

ORF Size: 1647 bp





SLC33A1 (NM_004733) Human Tagged Lenti ORF Clone - RC204131L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 004733.2</u>

RefSeq Size:3811 bpRefSeq ORF:1650 bpLocus ID:9197

 UniProt ID:
 000400

 Cytogenetics:
 3q25.31

Protein Families: Transmembrane

Protein Pathways: Glycosphingolipid biosynthesis - ganglio series, Metabolic pathways

MW: 60.9 kDa


Gene Summary: The protein encoded by this gene is required for the formation of O-acetylated (Ac)

gangliosides. The encoded protein is predicted to contain 6 to 10 transmembrane domains, and a leucine zipper motif in transmembrane domain III. Defects in this gene have been reported to cause spastic paraplegia autosomal dominant type 42 (SPG42) in one Chinese family, but not in similar patients of European descent. Two transcript variants encoding the


same protein have been found for this gene. [provided by RefSeq, Jul 2010]



Product images:



Circular map for RC204131L1



Double digestion of RC204131L1 using Sgfl and Mlul