

Product datasheet for **RC204133L1V**

STK38 (NM_007271) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | STK38 (NM_007271) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | STK38 |
| Synonyms: | NDR; NDR1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_007271 |
| ORF Size: | 1395 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC204133). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_007271.2 |
| RefSeq Size: | 3593 bp |
| RefSeq ORF: | 1398 bp |
| Locus ID: | 11329 |
| UniProt ID: | Q15208 |
| Cytogenetics: | 6p21.31 |
| Domains: | ptkase, TyrKc, S_TKc |
| Protein Families: | Druggable Genome, Protein Kinase |



[View online »](#)

MW: 54.2 kDa

Gene Summary: This gene encodes a member of the AGC serine/threonine kinase family of proteins. The kinase activity of this protein is regulated by autophosphorylation and phosphorylation by other upstream kinases. This protein has been shown to function in the cell cycle and apoptosis. This protein has also been found to regulate the protein stability and transcriptional activity of the MYC oncogene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2015]