

Product datasheet for **RC204146A1V**

Human FRA2 (FOSL2) (NM_005253) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FRA2 (FOSL2) (NM_005253) AAV Particle
Tag:	Myc-DDK
Symbol:	FRA2
Synonyms:	FRA2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204146 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTACCAGGATTATCCCGGAACTTTGACACCTCGTCCCGGGGAGCAGCGGCTCTCCTGCGCACGCCG
AGTCCTACTCCAGCGGCGGGCGGCCAGCAGAAATCCGGGTAGATATGCCTGGCTCAGGCAGTGCAATT
CATCCCCACCATCAACGCCATCACGACCAGCAGGACCTGCAGTGGATGGTGCAGCCACAGTGATCACC
TCCATGTCCAACCCATACCCCTCGCTCGCACCCCTACAGCCCCCTGCCGGGCTGGCCTCTGTCCCTGGAC
ACATGGCCCTCCCAAGACCTGGCGTGATCAAGACCATTGGCACCACCGTGGGCGCAGGAGGAGAGATGA
GCAGTGTCTCCTGAAGAGGAGGAGAAGCGTCGCATCCGGCGGGAGAGGAACAAGCTGGCTGCAGCCAAG
TGCCGGAACCGACGCGGGGAGCTGACAGAGAAGCTGCAGGCGGAGACAGAGGAGCTGGAGGAGGAGAAGT
CAGGCCTGCAGAAGGAGATTGCTGAGCTGCAGAAGGAGAAGGAGAAGCTGGAGTTTATGTTGGTGGCTCA
CGGCCAGTGTGCAAGATTAGCCCCGAGGAGCGCCGATCGCCCCAGCCCTGGGCTGCAGCCCATGCGC
AGTGGGGTGGCTCGGTGGGCGCTGTAGTGGTAAACAGGAGCCCTGGAAGAGGACAGCCCTCGTCCT
CGTCGGCGGGGCTGGACAAGGCCAGCGCTCTGTCAATCAAGCCATCAGCATTGCTGGGGCTTCTACGG
TGAGGAGCCCTGCACACCCCATCGTGGTGACCTCCACACCTGCTGCACTCCGGGCACCTCGAACCTC
GTCTTACCTATCCTAGCGTCTGGAGCAGGAGTACCCGCATCTCCCTCCGAATCCTGCTCCAAGGCTC
ACCGCAGAAGCAGTAGCAGCGGGACCAATCATCAGACTCCTTGAACCTCCCCACTCTGCTGGCTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204146 protein sequence Red=Cloning site Green=Tags(s) MYQDYPGNFDTSSRGSSGSPAHAESYSSGGGGQKFRVDMPGSGSAFIPTINAITTSQDLQWMVQPTVIT SMSNPYPRSHYPYSPGLASVPGHMLPRPGVIKTIGTTVGRRRRDEQLSPEEEKRRIRRRERKNLAAAK CRNRRRELTEKLQAETEELEEEKSGLQKEIAELQKEKEKLEFMLVAHGVPCKISPEERRSPAPGLQPMR SGGGSVGAVVVKQEPLEEDSPSSSSAGLDKAQRSVIKPIISAGGFYGEELHTPIVVTSTPAVTPGTSNL VFTYPSVLEQESPA SPSESCSAHRRSSSGDQSSDSLNSPTLLAL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005253
ORF Size:	978 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005253.3</u>
RefSeq Size:	4015 bp
RefSeq ORF:	981 bp
Locus ID:	2355
UniProt ID:	<u>P15408</u>
Cytogenetics:	2p23.2
MW:	35.2 kDa