

Product datasheet for RC204155L3

TM7SF1 (GPR137B) (NM_003272) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TM7SF1 (GPR137B) (NM_003272) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	TM7SF1
Synonyms:	TM7SF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204155).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_003272
ORF Size:	1197 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003272.1
RefSeq Size:	2046 bp
RefSeq ORF:	1200 bp
Locus ID:	7107
UniProt ID:	O60478
Cytogenetics:	1q42.3
Protein Families:	Druggable Genome, Transmembrane
MW:	45.6 kDa
Gene Summary:	Lysosomal integral membrane protein that regulates the localization and activity of mTORC1, a signaling complex promoting cell growth in response to growth factors, energy levels, and amino acids (PubMed:31036939). Interacts with Rag GTPases and increases the lysosomal localization and activity of Rag GTPases and thereby regulates mTORC1 translocation and activity in lysosome (PubMed:31036939). Involved in the regulation of lysosomal morphology and autophagy (PubMed:31036939).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC204155L3