

## Product datasheet for **RC204157L3V**

### **SNAP45 (SNAPC2) (NM\_003083) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SNAP45 (SNAPC2) (NM_003083) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNAP45
Synonyms:	PTFDELTA; SNAP45
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003083
ORF Size:	1002 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204157).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003083.1</a>
RefSeq Size:	1561 bp
RefSeq ORF:	1005 bp
Locus ID:	6618
UniProt ID:	<a href="#">Q13487</a>
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
MW:	35.6 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a subunit of the snRNA-activating protein complex which is associated with the TATA box-binding protein. The encoded protein is necessary for RNA polymerase II and III dependent small-nuclear RNA gene transcription. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2009]